



More information on the website
radwag.com/it/info,w1,DXV

Valigia di trasporto per Termobilancia MA X2, MAR



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri fisici

| | |
|----------------------|----------------|
| Dimensioni del pacco | 545×330×340 mm |
|----------------------|----------------|

Compatibile con

Termobilancia MA X7.IC.A
 Termobilancia MA 5Y
 Termobilancia MA X2.IC.A
 Termobilancia MA X2.A

Termobilancia MA X7.A
 Termobilancia MA 5Y.IC
 Termobilancia MA R

Device dimensions

