



More information on the website  
[radwag.com/en/info,w1,E8Z](http://radwag.com/en/info,w1,E8Z)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 500 g), Wooden Box

OK-301-0030



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	1 110 g
Nominal value	1 – 500 g
Quantity	12 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes