



More information on the website
radwag.com/en/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank



The drawings, photos and graphics used are for illustrative purposes only.

Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 WPY Multifunctional Scale
 Microbalance for Pipette Calibration MYA 5Y.P
 Tunnel Detectors
 C315 Platform scale
 5Y.PM Precision Balance
 Microbalance for Filters MYA 5Y.F
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balance
 C315.D Baby Scale
 XA 5Y.M.A Microbalance
 WLC Precision Balance
 PS R2 Precision Balance

AS 5Y Analytical Balance
 PS R2.H Precision Balance
 WTC Precision Balance
 CY10 Precision Balance
 C315.OW Personal Scale
 XA 5Y.F Analytical Balance for Filters
 WLY Precision Balance
 WLC C/2 Precision Balance
 MYA 5Y Microbalance
 PS R1 Precision Balance
 XA 5Y Analytical Balance
 AS R2 PLUS Analytical Balance
 PUE C32 Weighing Terminal
 PS 5Y Precision Balance

WLC X7 Precision Balance
PUE CY10 Weighing Terminal
C315.4 Platform Scale
C32.PM Precision Balance
C315.0 Personal Scale
PS X7 Precision Balance
XA 5Y.A Analytical Balance
C315.4P pallet scales
UYA 5Y Ultra-Microbalance
Powerbanks

XA 5Y.M Microbalance
PUE C315 Weighing Terminal
AS R1 PLUS Analytical Balance
Ultra-Microbalance for Filters UYA 5Y.F
C315.4B Bed Scale
C315.K Chair Scale
XA 5Y.M.A.S Microbalance for stents weighing
C315.8B Ramps Bed Scale
Receipt Printer