



More information on the website  
[radwag.com/en/info,w1,F64](https://radwag.com/en/info,w1,F64)

## Cable to the ZR-03/1 and ZR-04/1 Powerbank



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

### Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
WPY Multifunctional Scale  
MYA 5Y.P Microbalance for pipettes calibration  
Tunnel Detectors  
C315 Platform scale  
5Y.PM Precision Balance  
MYA 5Y.F Microbalance for Filters  
AS X7 Analytical Balance  
PS X2 Precision Balance  
AS X2 PLUS Analytical Balance  
C315.D Baby Scale  
XA 5Y.M.A Microbalance  
WLC Precision Balance  
PS R2 Precision Balance

PS 5Y Balances for PGC  
PS R2.H Precision Balance  
WTC Precision Balance  
CY10 Precision Balance  
C315.OW Personal Scale  
XA 5Y.F Analytical Balance for Filters  
WLY Precision Balance  
WLC C/2 Precision Balance  
MYA 5Y Microbalance  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
AS R2 PLUS Analytical Balance  
PUE C32 Weighing Terminal  
PS 5Y Precision Balance

WLC X7 Precision Balance  
PUE CY10 Weighing Terminal  
C315.4 Platform Scale  
C32.PM Precision Balance  
C315.0 Personal Scale  
PS X7 Precision Balance  
XA 5Y.A Analytical Balance  
C315.4P pallet scales  
UYA 5Y Ultra-Microbalance  
AS 5Y Analytical Balance

XA 5Y.M Microbalance  
PUE C315 Weighing Terminal  
AS R1 PLUS Analytical Balance  
UYA 5Y.F Ultra-Microbalance for Filters  
C315.4B Bed Scale  
C315.K Chair Scale  
XA 5Y.M.A.S Microbalance for stents weighing  
C315.8B Ramps Bed Scale  
Receipt Printer