



More information on the website
radwag.com/en/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank

WX-013-0244



The drawings, photos and graphics used are for illustrative purposes only.

Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
WPY Multifunctional Scale
MYA 5Y.P Microbalance for pipettes calibration
Tunnel Detectors
C315 Platform scale
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balance
C315.D Baby Scale
XA 5Y.M.A Microbalance
WLC Precision Balance

PS 5Y Balances for PGC
PS R2.H Precision Balance
WTC Precision Balance
CY10 Precision Balance
C315.OW Personal Scale
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
WLC C/2 Precision Balance
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal

PS R2 Precision Balance
WLC X7 Precision Balance
PUE CY10 Weighing Terminal
C315.4 Platform Scale
C32.PM Precision Balance
C315.0 Personal Scale
PS X7 Precision Balance
XA 5Y.A Analytical Balance
C315.4P pallet scales
UYA 5Y Ultra-Microbalance
AS 5Y Analytical Balance

PS 5Y Precision Balance
XA 5Y.M Microbalance
PUE C315 Weighing Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed Scale
C315.K Chair Scale
XA 5Y.M.A.S Microbalance for stents weighing
C315.8B Ramps Bed Scale
Receipt Printer