



More information on the website
radwag.com/us/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank



The drawings, photos and graphics used are for illustrative purposes only.

Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatible with

XA 5Y.M.A.P Microbalance
WPY Multifunctional Scale
Microbalance for Pipette Calibration MYA 5Y.P
Tunnel detectors
C315 Platform scale
5Y.PM Precision balance
Microbalance for Filters MYA 5Y.F
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balance
C315.D Baby scale
XA 5Y.M.A Microbalance
WLC Precision Balance
PS R2 Precision Balance

AS 5Y Analytical Balance
PS R2.H Precision Balance
WTC Precision Balance
CY10 Precision balance
C315.OW Personal scale
Analytical Balance for Filters XA 5Y.F
WLY Precision Balance
WLC C/2 Precision balance
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PS 5Y Precision Balance

WLC X7 Precision balance
PUE CY10 Weighing Terminal
C315.4 Platform Scale
C32.PM Precision Balance
C315.O Personal scale
PS X7 Precision Balance
XA 5Y.A Analytical Balance
C315.4P pallet scales
UYA 5Y Ultra-Microbalance
Powerbanks

XA 5Y.M Microbalance
PUE C315 Terminal
AS R1 PLUS Analytical Balance
Ultra-Microbalance for Filters UYA 5Y.F
C315.4B Bed scale
C315.K Chair scale
XA 5Y.M.A.S Microbalance
C315.8B Ramps bed scale
Receipt Printer