



More information on the website
radwag.com/us/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank

WX-013-0244



The drawings, photos and graphics used are for illustrative purposes only.

Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 WPY Multifunctional Scale
 MYA 5Y.P Microbalance for pipettes calibration
 Tunnel detectors
 C315 Platform scale
 5Y.PM Precision balance
 MYA 5Y.F Microbalance for Filters
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balance
 C315.D Baby scale
 XA 5Y.M.A Microbalance
 WLC Precision Balance

PS 5Y Balances for PGC
 PS R2.H Precision Balance
 WTC Precision Balance
 CY10 Precision balance
 C315.OW Personal scale
 Analytical Balance for Filters XA 5Y.F
 WLY Precision Balance
 WLC C/2 Precision balance
 MYA 5Y Microbalance
 PS R1 Precision Balance
 XA 5Y Analytical Balance
 AS R2 PLUS Analytical Balance
 PUE C32 Weighing Terminal

PS R2 Precision Balance
WLC X7 Precision balance
PUE CY10 Weighing Terminal
C315.4 Platform Scale
C32.PM Precision Balance
C315.0 Personal scale
PS X7 Precision Balance
XA 5Y.A Analytical Balance
C315.4P pallet scales
UYA 5Y Ultra-Microbalance
AS 5Y Analytical Balance

PS 5Y Precision Balance
XA 5Y.M Microbalance
PUE C315 Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed scale
C315.K Chair scale
XA 5Y.M.A.S Microbalance for stents weighing
C315.8B Ramps bed scale
Receipt Printer