



More information on the website
radwag.com/en/info,w1,FBE

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 2 kg), Plastic Box

OK-202-0069



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	6 111.11 g
Nominal value	1 mg – 2 kg
Quantity	27 pcs

Compatible with (Additional Fee)

Tweezers

Brushes