



More information on the website
radwag.com/en/info,w1,FOT

USB 3.0 Hub

WX-015-0286



The drawings, photos and graphics used are for illustrative purposes only.

Overview

The HUB is powered via the USB and does not require driver installation. Simply connect the USB cable to your device and computer and start working. The HUB features surge protection in the event of a short circuit on peripheral devices.

The HUB is equipped with 4 USB 3.0 Plug & Play ports. They make it possible to connect enough necessary devices. The USB 3.0 ports are particularly important for optimum data transfer from a USB 3.0 flash drive.



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
WPY Multifunctional Scale
MYA 5Y.P Microbalance for pipettes calibration
MA 5Y Moisture Analyzer
5Y.PM Precision Balance

XA 5Y.A Analytical Balance
MA 5Y.IC Moisture Analyzer
AS 5Y Analytical Balance
PS 5Y Balances for PGC
CY10 Precision Balance

MYA 5Y.F Microbalance for Filters
WLY Scales for PGC
XA 5Y.M.A Microbalance
CY10 Scales for PGC
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
PUE 7.1.HRP High Resolution Scale
Automatic Balances for. Pipettes Calibration
WPY Scales for PGC

XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
MYA 5Y Microbalance
XA 5Y Analytical Balance
PUE 7.1 Weighing Terminal
PS 5Y Precision Balance
XA 5Y.M Microbalance
H315 H Waterproof Platform Scale
XA 5Y.M.A.S Microbalance for stents weighing