



More information on the website  
[radwag.com/it/info,w1,FWC](http://radwag.com/it/info,w1,FWC)

# Pesi standard F1 - cilindrico - 500 g

OP-203-0120



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

| Parametri metrologici |                 |
|-----------------------|-----------------|
| Classe OIML           | F1              |
| Costruzione           |                 |
| Materiale             | Stainless steel |
| Valore nominale       | 500 g           |

## Compatibile con (Additional Fee)

Pennelli

Pinzette