



# Pesi standard F2 - cilindrico - 200 g

OP-303-0053

More information on the website  
[radwag.com/it/info,w1,GR5](http://radwag.com/it/info,w1,GR5)



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

|             |    |
|-------------|----|
| Classe OIML | F2 |
|-------------|----|

### Costruzione

|           |                 |
|-----------|-----------------|
| Materiale | Stainless steel |
|-----------|-----------------|

|                 |       |
|-----------------|-------|
| Valore nominale | 200 g |
|-----------------|-------|

## Compatibile con (Additional Fee)

Pinzette

Pennelli