



More information on the website  
[radwag.com/en/info,w1,GZV](https://radwag.com/en/info,w1,GZV)

# SAL/H/PLUS Stainless Steel Laboratory Antivibration Table for PLUS Series Balances



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	77 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 MYA 5Y.F Microbalance for Filters

MYA 5Y Microbalance  
 XA 5Y Analytical Balance  
 XA 5Y.KO Manual Mass Comparator

XA 5Y.M.A Microbalance  
XA 5Y.A Analytical Balance  
UYA 5Y Ultra-Microbalance  
AS 5Y Analytical Balance  
XA 5Y.F Analytical Balance for Filters

XA 5Y.M Microbalance  
UYA 5Y.KO Manual Mass Comparator  
UYA 5Y.F Ultra-Microbalance for Filters  
XA 5Y.M.A.S Microbalance for stents weighing

## Device dimensions

SAL/H/PLUS

