



More information on the website
radwag.com/en/info,w1,GZV

SAL/H/PLUS Stainless Steel Laboratory Antivibration Table for PLUS Series Balances



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	77 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters

XA 5Y.F Analytical Balance for Filters
 MYA 5Y Microbalance
 XA 5Y Analytical Balance

XA 5Y.M.A Microbalance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance

XA 5Y.KO Manual Mass Comparator
XA 5Y.M Microbalance
XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions

SAL/H/PLUS

