



More information on the website
radwag.com/en/info,w1,GZV

SAL/H/PLUS Stainless Steel Laboratory Antivibration Table for PLUS Series Balances

WX-001-0155



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	77 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration

MYA 5Y Microbalance

MYA 5Y.P Microbalance for pipettes calibration
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
PS R2 Precision Balance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
PS R2.H Precision Balance
XA 5Y.F Analytical Balance for Filters

PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
XA 5Y.KO Manual Mass Comparator
XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H

SAL/H/PLUS

