



More information on the website
radwag.com/us/info,w1,GZV

SAL/H/PLUS – Anti-vibration table in stainless steel technology



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	77 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters
 XA 5Y.M.A Microbalance

MYA 5Y Microbalance
 XA 5Y Analytical Balance
 XA 5Y.KO Manual Mass Comparator
 XA 5Y.M Microbalance

XA 5Y.A Analytical Balance
UYA 5Y Ultra-Microbalance
AS 5Y Analytical Balance
Analytical Balance for Filters XA 5Y.F

UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions

SAL/H/PLUS

