



More information on the website  
[radwag.com/de/info,w1,HCJ](http://radwag.com/de/info,w1,HCJ)

# Prüfgewicht der M1 Klasse - Knopfgewicht - 1 g

OP-403-0118



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Konstruktion	
Material	Edelstahl
Nominalwert	1 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself