



More information on the website
radwag.com/us/info,w1,HHO

PT0454 Cable



The drawings, photos and graphics used are for illustrative purposes only.

Overview

RS 232 cable for connecting a balance/scale to a computer.

Cable length: 5 m

Compatible with:

- Scales operated using the PUE C315 terminal
- PUE C315 terminal

Compatible with

XA 5Y.M.A.P Microbalance

MediaBox

Microbalance for Pipette Calibration MYA 5Y.P

5Y.PM Precision balance

Microbalance for Filters MYA 5Y.F

XA 5Y.M.A Microbalance

PUE CY10 Weighing Terminal

CY10 Multifunctional Scale

XA 5Y.A Analytical Balance

UYA 5Y Ultra-Microbalance

CY10 Precision balance

Analytical Balance for Filters XA 5Y.F

MYA 5Y Microbalance

XA 5Y Analytical Balance

PMV 5Y Moisture Analyzer

XA 5Y.M Microbalance

Ultra-Microbalance for Filters UYA 5Y.F

XA 5Y.M.A.S Microbalance