



More information on the website
radwag.com/en/info,w1,IFA

Knob weights, Set (1 mg - 200 g), Wooden Box

OK-101-0042



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | E2 |
| Construction | |
| Material | Stainless steel |
| Weight | 611.11 g |
| Nominal value | 1 mg – 200 g |
| Quantity | 23 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes