



More information on the website  
[radwag.com/en/info,w1,KAA](https://radwag.com/en/info,w1,KAA)

# Knob weights, Set (1 mg - 2 kg), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 6 111.11 g      |
| Nominal value           | 1 mg – 2 kg     |
| Quantity                | 27 pcs          |

## Compatible with

Tweezers

Brushes