



More information on the website  
[radwag.com/en/info,w1,KD6](http://radwag.com/en/info,w1,KD6)

# Tweezers 250 mm for mass standards up to 500 g



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Material	Stainless steel
Physical parameters	
Length	250 mm

## Compatible with

E1 Mass Standard - knob weights, sets (1 mg - 10 kg)  
 F1 Mass Standard - stainless steel wires  
 F1 Mass Standard - knob weights with adjustment chamber, sets  
 F1 Mass Standard - knob weights without adjustment chamber  
 F2 Mass Standard - stainless steel wires  
 F2 Mass Standard - Carrier weights  
 E2 Mass Standard - stainless steel wires  
 E2 Mass Standard - knob weights, sets (1 mg - 10 kg)  
 F2 Mass Standard - knob weights with adjustment chamber, sets  
 E2 Mass Standard - cylindrical weights, sets (1 g - 10 kg)  
 F2 Mass Standard - Slotted weights  
 M1 Mass Standard - flat polygonal sheets  
 F2 Mass Standard - flat polygonal sheets  
 F1 Mass Standard - knob weights without adjustment chamber, set

E2 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)  
 F1 Mass Standard - knob weights with adjustment chamber  
 E2 Mass Standard - knob weights  
 E2 Mass Standard - flat polygonal sheets  
 M1 Mass Standard - knob weights  
 F1 Mass Standard - flat polygonal sheets, sets  
 E1 Mass Standard - knob weights  
 E1 Mass Standard - stainless steel wires  
 F2 Mass Standard - knob weights with adjustment chamber  
 F2 Mass Standard - flat polygonal sheets, sets  
 E1 Mass Standard - flat polygonal sheets  
 E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)  
 E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)  
 F2 Mass Standard - Stainless steel wires, sets

F1 Mass Standard - flat polygonal sheets  
E1 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)  
F1 Mass Standard - cylindrical weights

F1 Mass Standard - stainless steel wires, sets  
M1 Mass Standard - flat polygonal sheets, sets