



More information on the website
radwag.com/en/info,w1,KD6

Tweezers 250 mm for mass standards up to 500 g



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Material	Stainless steel
Physical parameters	
Length	250 mm

Compatible with

E1 Mass Standard - knob weights, sets (1 mg - 10 kg)
 F1 Mass Standard - stainless steel wires
 F1 Mass Standard - knob weights with adjustment chamber, sets
 F1 Mass Standard - knob weights without adjustment chamber
 F2 Mass Standard - stainless steel wires
 F2 Mass Standard - Carrier weights
 E2 Mass Standard - stainless steel wires
 E2 Mass Standard - knob weights, sets (1 mg - 10 kg)
 F2 Mass Standard - knob weights with adjustment chamber, sets
 E2 Mass Standard - cylindrical weights, sets (1 g - 10 kg)
 F2 Mass Standard - Slotted weights
 M1 Mass Standard - flat polygonal sheets
 F2 Mass Standard - flat polygonal sheets
 F1 Mass Standard - knob weights without adjustment chamber, set

E2 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
 F1 Mass Standard - knob weights with adjustment chamber
 E2 Mass Standard - knob weights
 E2 Mass Standard - flat polygonal sheets
 M1 Mass Standard - knob weights
 F1 Mass Standard - flat polygonal sheets, sets
 E1 Mass Standard - knob weights
 E1 Mass Standard - stainless steel wires
 F2 Mass Standard - knob weights with adjustment chamber
 F2 Mass Standard - flat polygonal sheets, sets
 E1 Mass Standard - flat polygonal sheets
 E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
 E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
 F2 Mass Standard - Stainless steel wires, sets

F1 Mass Standard - flat polygonal sheets
E1 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
F1 Mass Standard - cylindrical weights

F1 Mass Standard - stainless steel wires, sets
M1 Mass Standard - flat polygonal sheets, sets