



More information on the website  
[radwag.com/en/info,w1,KUV](http://radwag.com/en/info,w1,KUV)

# F2 Mass Standard - Knob Weights With Adjustment Chamber - 500 g



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F2              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 500 g           |

## Compatible with

Brushes

Tweezers