



More information on the website  
[radwag.com/en/info,w1,LOG](https://radwag.com/en/info,w1,LOG)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 500 g), Plastic Box

OK-302-0037



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	1 111.11 g
Nominal value	1 mg - 500 g
Quantity	24 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes