



More information on the website
radwag.com/en/info,w1,LXB

DS2208 Barcode Scanner



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Protection class	IP 42
Communication interface	
Communication interface	RS232, USB
Environmental conditions	
Operating temperature	0 ÷ +50 °C
Storage temperature	-40 ÷ +70 °C
Physical parameters	
Device dimensions	165×99×66 mm
Drop resistance [m]:	1.5
Operated barcodes:	1D, 2D, Pdf417, Postal, GS1 Databar
Readout technology:	Imager 2D (area imager)
Required code contrast [%]:	25

Permissible ambient humidity [%]:

5% ÷ 95%

Compatible with

MUYA Modules

HX7 H Stainless steel Multifunctional Scale
HX7.4N H Multifunctional Stainless Steel Ramp Scale
HX7.4C Multifunctional Platform Scale
HX7 HR Multifunctional Scale
XA 5Y.M.A.P Microbalance for pipettes calibration
WPY Multifunctional Scale
Microbalance for Pipette Calibration MYA 5Y.P
C32 Multifunctional Scale
HY10.HRP.H Stainless steel High Resolution Scale
5Y.PM Precision Balance
Microbalance for Filters MYA 5Y.F
HX7.4P2 H Multifunctional Stainless Steel Beam Scales
TMX H Stainless steel Multifunctional Scale
WLY Scales for PGC
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balance
TMX Multifunctional Scale
XA 5Y.M.A Microbalance
PS R2 Precision Balance
WLC X7 Precision Balance
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale
PS X2 Balances for PGC
DWT/RC/HYF Checkweighers
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
HX7.4P H Multifunctional Stainless Steel Pallet Scales
HY10 Multifunctional Scale
HY10 H Stainless Steel Multifunctional Scale
HY10.HRP High Resolution Scale
C32.PM Precision Balance
PUE 7.1.HRP High Resolution Scale

Automatic Balances for. Pipettes Calibration

HX7 Multifunctional Scale
HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version
PS X7 Precision Balance
WPY Scales for PGC
XA 5Y.A Analytical Balance
HX7.4 H Multifunctional Stainless Steel Platform Scale
HX7.4P2 Multifunctional Beam Scales
WAY 5Y.KO Manual Mass Comparator
AS 5Y Analytical Balance
PS R2.H Precision Balance
HY10.H Scales for PGC
CY10 Precision Balance
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
APP 5Y.KO Manual Mass Comparator
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PUE HX7 Weighing Terminal
PUE 5 Weighing Terminal
PS 5Y Precision Balance
XA 5Y.KO Manual Mass Comparator
XA 5Y.M Microbalance
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
H315 H Waterproof Platform Scale
HX7.4P Multifunctional pallet scales
XA 5Y.M.A.S Microbalance for stents weighing