



More information on the website  
[radwag.com/en/info,w1,M9A](http://radwag.com/en/info,w1,M9A)

# Knob weights, Set (1 kg - 5 kg), Wooden Box

OK-101-0038



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E2              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 10 000 g        |
| Nominal value           | 1 – 5 kg        |
| Quantity                | 4 pcs           |

## Compatible with (Additional Fee)

Brushes