



More information on the website  
[radwag.com/en/info,w1,MTS](https://radwag.com/en/info,w1,MTS)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 1 kg), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	2 110 g
Nominal value	1 g - 1 kg
Quantity	13 pcs

## Compatible with

Tweezers

Brushes