



More information on the website
radwag.com/it/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances



The drawings, photos and graphics used are for illustrative purposes only.

Foglio di calcolo

Parametri fisici	
Dimensioni	560x350x352 mm
Dimensioni del pacco	695x485x485 mm
Peso netto	11 kg
Peso lordo	16 kg

Compatibile con

Microbalance XA 5Y.M.A.P
 Microbalance XA 5Y.M.A
 Bilance analitiche XA 5Y.A
 Bilance analitiche XA 5Y


XA 5Y.KO Comparatore di massa manuale
 Microbalance XA 5Y.M
 Microbalance XA 5Y.M.A.S

Device dimensions




Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances



Anti-draft chamber for microbalances




Draft shield for precision balances



Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances