



More information on the website
radwag.com/us/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet


Physical parameters	
Dimension	560x350x352 mm
Packaging dimensions	695x485x485 mm
Net weight	11 kg
Gross weight	16 kg

Compatible with

XA 5Y.M.A.P Microbalance
 XA 5Y.M.A Microbalance
 XA 5Y.A Analytical Balance
 XA 5Y Analytical Balance


XA 5Y.KO Manual Mass Comparator
 XA 5Y.M Microbalance
 XA 5Y.M.A.S Microbalance

Device dimensions




Anti-draft chamber for balances with 128x128 mm weighing pan


Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances