



More information on the website  
[radwag.com/us/info,w1,ODT](http://radwag.com/us/info,w1,ODT)

# DJ-04 Antistatic ionizer



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Power supply	100 – 240 V AC 50/60 Hz
Operating temperature	0 ÷ +40 °C
Relative humidity	5% ÷ 80%
Distance form a sample	Ok. 5 ÷ 50 mm
Ion source	stainless steel
<b>Physical parameters</b>	
Frame dimensions	268×110×295 mm
Packaging dimensions	375×375×130 mm
Net weight	1,3 kg
Gross weight	2,5 kg

## Compatible with


XA 5Y.M.A.P Microbalance  
 AS X7 Analytical Balance

Analytical Balance for Filters XA 5Y.F  
 XA 5Y Analytical Balance

AS X2 PLUS Analytical Balances  
XA 5Y.M.A Microbalance  
XA 5Y.A Analytical Balance  
AS 5Y Analytical Balance


AS R2 PLUS Analytical Balance  
XA 5Y.M Microbalance  
AS R1 PLUS Analytical Balance  
XA 5Y.M.A.S Microbalance

## Device dimensions




DJ-04

---



DJ-05

---



power supply

---