



More information on the website  
[radwag.com/en/info,w1,P1D](http://radwag.com/en/info,w1,P1D)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 1 kg), Plastic Box

OK-202-0076



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	2 110 g
Nominal value	1 g – 1 kg
Quantity	13 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes