



More information on the website
radwag.com/de/info,w1,PZ0

Prüfgewicht der M1 Klasse - Knopfgewicht - 10 g

OP-403-0121



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | M1 |
| Konstruktion | |
| Material | Edelstahl |
| Nominalwert | 10 g |

Kompatibel mit (Additional Fee)

Pinzetten

Pinzel