



More information on the website  
radwag.com/us/info,w1,QKR

# Stand for terminal 5Y / PUE CY10 - mild steel version

WX-004-0298



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A stand, based on a tripod for supporting scale's indicator installed on a cable.

## Compatible with (Additional Fee)

5Y.PM Precision balance  
CY10 Scales for PGC  
CY10 Multifunctional Scale

CY10 Precision balance  
HRP 5Y.KO Mass Comparator  
HRP 5Y.KB Mass Comparator

## Device dimensions W x D x H

