



More information on the website
radwag.com/us/info,w1,QKR

Stand for terminal 5Y / PUE CY10 - mild steel version



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A stand, based on a tripod for supporting scale's indicator installed on a cable.

Compatible with

5Y.PM Precision balance
CY10 Scales for PGC
CY10 Multifunctional Scale

CY10 Precision balance
HRP 5Y.KO Mass Comparator
HRP 5Y.KB Mass Comparator

Device dimensions

