



More information on the website  
[radwag.com/en/info,w1,QMN](http://radwag.com/en/info,w1,QMN)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 1 kg), Plastic Box

OK-302-0038



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	2 111.11 g
Nominal value	1 mg – 1 kg
Quantity	25 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes