



More information on the website
radwag.com/us/info,w1,QNI

Voluminous sample container

WX-003-0029



The drawings, photos and graphics used are for illustrative purposes only.

Overview

The set is an ideal solution for drying large volume samples, e.g. textiles, or products that placed on the weighing pan tend to increase their volume. The container allows to enclose the sample within a limited space wherein the possibility to determine the moisture content value remains unaffected.

In the case of large-volume samples, the kit prevents the sample from coming into contact with solid elements of the moisture analyzer (such as the heating element or its cover) and protects it from damage. Closing a large-volume, low-density sample in the container also reduces the risk of errors resulting from the effect of the sample being raised by warm air.

The accessory gives the user full control of the sample and an error-free result.

The set consists of 2 components made of steel. It is compatible with all moisture analyzers.

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MA 5Y Moisture Analyzer
MA X7.A Moisture Analyzer

MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer
MA X2.A Moisture Analyzers