



More information on the website
radwag.com/en/info,w1,RG1

F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 50 g), Plastic Box



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 111.11 g |
| Nominal value | 1 mg – 50 g |
| Quantity | 20 pcs |

Compatible with

WPY Multifunctional Scale
 TMX H Stainless steel Multifunctional Scale
 WLY Scales for PGC
 HY10 Multifunctional Scale
 PUE 5 Weighing Terminal
 HY10 H Stainless Steel Multifunctional Scale

WPY Scales for PGC
 HY10.H Scales for PGC
 Tweezers
 WLY Precision Balance
 Brushes