



More information on the website  
[radwag.com/en/info,w1,RTH](https://radwag.com/en/info,w1,RTH)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 1 kg), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F2              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 2 110 g         |
| Nominal value           | 1 g - 1 kg      |
| Quantity                | 13 pcs          |

## Compatible with

Tweezers

Brushes