



More information on the website  
[radwag.com/de/info,w1,S20](http://radwag.com/de/info,w1,S20)

# Knopfgewichte, Set (1 g - 500 g), Kunststoff-Etui

OK-402-0067



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Messtechnische Parameter

OIML-Klasse	M1
Gewicht	1 110 g
Nominalwert	1 – 500 g
Anzahl	12 St.

## Kompatibel mit (Additional Fee)

Pinsel