



More information on the website
radwag.com/en/info,w1,SCH

P0108 Cable



The drawings, photos and graphics used are for illustrative purposes only.

Overview

RS 232 cable: Balance - computer. Standard cable length: 1.7 m. Compatible with:

- Scales with PUE C315 terminal
- PUE C315 terminal

Compatible with

MA X7.IC.A Moisture Analyzer
 MUYA Modules
 WPY Multifunctional Scale
 C32 Multifunctional Scale
 MA 5Y Moisture Analyzer
 C315 Platform scale
 WLY Scales for PGC
 MA X2.IC.A Moisture Analyzer
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balance
 PS X7 Balances for PGC
 C315.D Baby Scale
 WLC Precision Balance

PS 5Y Balances for PGC
 AKM-2 Automatic Mass Comparator
 PS R2.H Precision Balance
 WTC Precision Balance
 H315 4C Platform Scale
 C315.OW Personal Scale
 XA 5Y.F Analytical Balance for Filters
 WLY Precision Balance
 UMA Automatic Mass Comparator
 APP 5Y.KO Manual Mass Comparator
 Robotic Mass Comparators
 Density Measurement
 WLC C/2 Precision Balance
 PS R1 Precision Balance

PS R2 Precision Balance
C315.4P2 Beam Scales
WLC X7 Precision Balance
PS X2 Balances for PGC
MA X2.A Moisture Analyzer
MAS Modules
MPS Modules
C315.4 Platform Scale
C32.PM Precision Balance
MA X7.A Moisture Analyzer
C315.0 Personal Scale
PS X7 Precision Balance
WPY Scales for PGC
C315.4P pallet scales
WAY 5Y.KO Manual Mass Comparator
MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator

MA X2.A Moisture Analyzers
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PS 5Y Precision Balance
Postal Scales for Packages
Automatic Vacuum Mass Comparators
MAS1 Modules
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
PUE C315 Weighing Terminal
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
C315.4B Bed Scale
Nano-Comparator NANO.AK-4.500.5Y
C315.K Chair Scale
C315.8B Ramps Bed Scale