



More information on the website
radwag.com/en/info,w1,U1A

SAL/H Stainless Steel Laboratory Antivibration Table



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	65-70 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with

MUYA Modules
 MAS1 Modules
 MAS Modules

MPS Modules
 C32.PM Precision Balance
 XA 5Y.F Analytical Balance for Filters

Device dimensions

SAL/H

