



More information on the website
radwag.com/en/info,w1,U1A

SAL/H Stainless Steel Laboratory Antivibration Table

WX-001-0077



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	65-70 kg
Length	770 mm
Width	680 mm
Height	815±10 mm



Compatible with (Additional Fee)

AS X7 Analytical Balance
PS R2 Precision Balance
WLC X7 Precision Balance
WLC Precision Balance
C32.PM Precision Balance
PS R2.H Precision Balance
WTC Precision Balance

XA 5Y.F Analytical Balance for Filters
PS X7 Precision Balance
WLC C/2 Precision Balance
PS R1 Precision Balance
AS R2 PLUS Analytical Balance
AS R1 PLUS Analytical Balance
MAS.1.Y Modules

Device dimensions W x D x H

SAL/H

