



More information on the website
radwag.com/en/info,w1,U3R

PM9300-DK Barcode Scanner

WX-006-0034



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Protection class	IP 65
Environmental conditions	
Operating temperature	-20 – +50 °C
Storage temperature	-40 – +70 °C
Physical parameters	
Device dimensions	212×110×74 mm
Drop resistance [m]:	2
Operated barcodes:	1D, GS1 Databar
Readout technology:	Single line laser
Required code contrast [%]:	25
Metrological parameters	
Permissible ambient humidity [%]:	0 % – 95 %

Compatible with (Additional Fee)

MUYA Modules

HX7 H Stainless steel Multifunctional Scale
HX7.4N H Multifunctional Stainless Steel Ramp Scale
HX7.4C Multifunctional Platform Scale
HX7 HR Multifunctional Scale
XA 5Y.M.A.P Microbalance for pipettes calibration
WPY Multifunctional Scale
MYA 5Y.P Microbalance for pipettes calibration
C32 Multifunctional Scale
HY10.HRP.H Stainless steel High Resolution Scale
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
HX7.4P2 H Multifunctional Stainless Steel Beam Scales
TMX H Stainless steel Multifunctional Scale
WLY Scales for PGC
AS X7 Analytical Balance
PS X7 Balances for PGC
TMX Multifunctional Scale
XA 5Y.M.A Microbalance
PS R2 Precision Balance
WLC X7 Precision Balance
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale
CY10 Scales for PGC
DWT/RC/HYF Checkweighers
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
HX7.4P H Multifunctional Stainless Steel Pallet Scales
HY10 Multifunctional Scale
HY10 H Stainless Steel Multifunctional Scale
HY10.HRP High Resolution Scale
C32.PM Precision Balance
PUE 7.1.HRP High Resolution Scale
Automatic Balances for. Pipettes Calibration
HX7 Multifunctional Scale

HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version
PS X7 Precision Balance
WPY Scales for PGC
XA 5Y.A Analytical Balance
HX7.4 H Multifunctional Stainless Steel Platform Scale
HX7.4P2 Multifunctional Beam Scales
WAY 5Y.KO Manual Mass Comparator
AS 5Y Analytical Balance
PS 5Y Balances for PGC
PS R2.H Precision Balance
HY10.H Scales for PGC
CY10 Precision Balance
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
APP 5Y.KO Manual Mass Comparator
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PUE HX7 Weighing Terminal
PUE 5 Weighing Terminal
PS 5Y Precision Balance
XA 5Y.KO Manual Mass Comparator
XA 5Y.M Microbalance
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
H315 H Waterproof Platform Scale
HX7.4P Multifunctional pallet scales
DWM Checkweighers
Roller Scales
XA 5Y.M.A.S Microbalance for stents weighing