



More information on the website  
radwag.com/us/info,w1,U4E

# Suction pump for liquids

WX-001-0445



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

Precision suction, error-free measurements The product is supplied with two internal pumps and two built-in buttons used to operate separate pumps. Backlit buttons inform you about operation of a specific pump. The set has two tips for sucking water out of a single beaker (one-channel tip) or four pipettes at a time (four-channel tip). Water is carried from both channels to the beaker through one tube. The power is provided to the pump through a USB A in the weighing terminal. You can supply power from the phone charger, for instance.

## Compatible with (Additional Fee)

Automatic Balances for. Pipettes Calibration  
XA 5Y.M.A Microbalance

XA 5Y Analytical Balance  
XA 5Y.M Microbalance

