



More information on the website
radwag.com/en/info,w1,UYE

Tweezers 300 mm for mass standards up to 500 g



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Construction | |
|---------------------|-----------------|
| Material | Stainless steel |
| Physical parameters | |
| Length | 300 mm |

Compatible with

- E1 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F1 Mass Standard - stainless steel wires
- F1 Mass Standard - knob weights with adjustment chamber, sets
- F1 Mass Standard - knob weights without adjustment chamber
- F2 Mass Standard - stainless steel wires
- F2 Mass Standard - Carrier weights
- E2 Mass Standard - stainless steel wires
- E2 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F2 Mass Standard - knob weights with adjustment chamber, sets
- E2 Mass Standard - cylindrical weights, sets (1 g - 10 kg)
- F2 Mass Standard - Slotted weights
- M1 Mass Standard - flat polygonal sheets
- F2 Mass Standard - flat polygonal sheets
- F1 Mass Standard - knob weights without adjustment chamber, set

- E2 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
- F1 Mass Standard - knob weights with adjustment chamber
- E2 Mass Standard - knob weights
- E2 Mass Standard - flat polygonal sheets
- M1 Mass Standard - knob weights
- F1 Mass Standard - flat polygonal sheets, sets
- E1 Mass Standard - knob weights
- E1 Mass Standard - stainless steel wires
- F2 Mass Standard - knob weights with adjustment chamber
- F2 Mass Standard - flat polygonal sheets, sets
- E1 Mass Standard - flat polygonal sheets
- E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- F2 Mass Standard - Stainless steel wires, sets

F1 Mass Standard - flat polygonal sheets
E1 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
F1 Mass Standard - cylindrical weights

F1 Mass Standard - stainless steel wires, sets
M1 Mass Standard - flat polygonal sheets, sets