



More information on the website  
[radwag.com/en/info,w1,VDB](https://radwag.com/en/info,w1,VDB)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber - 200 g



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 200 g           |

## Compatible with

Tweezers

Brushes