



More information on the website  
[radwag.com/it/info,w1,VDB](http://radwag.com/it/info,w1,VDB)

## Pesi standard F1 - cilindrico - 200 g



The drawings, photos and graphics used are for illustrative purposes only.

### Foglio di calcolo

| Parametri metrologici |                 |
|-----------------------|-----------------|
| Classe OIML           | F1              |
| Costruzione           |                 |
| Materiale             | Stainless steel |
| Valore nominale       | 200 g           |

### Compatibile con

Pinzette

Pennelli