



More information on the website  
[radwag.com/en/info,w1,VIW](http://radwag.com/en/info,w1,VIW)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 2 kg), Plastic Box

OK-202-0077



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	6 110 g
Nominal value	1 g – 2 kg
Quantity	15 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes