



More information on the website
radwag.com/en/info,w1,W20

Under-pan Weighing Rack for AS and PS Balances

WX-002-0111



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-pan weighing intended for AS and PS balances.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with (Additional Fee)

AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balances
PS X7 Balances for PGC
PS R2 Precision Balance
PS X2 Balances for PGC
WLC Precision Balance
PS X7 Precision Balance

AS 5Y Analytical Balance
PS 5Y Balances for PGC
PS R2.H Precision Balance
WLC C/2 Precision Balance
PS R1 Precision Balance
AS R2 PLUS Analytical Balance
PS 5Y Precision Balance
AS R1 PLUS Analytical Balance

Device dimensions

