



More information on the website
radwag.com/us/info,w1,W20

Under-hook weighing rack for AS and PS balances

WX-002-0111



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-hook weighing intended for balances AS and PS.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with (Additional Fee)

AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balances
 PS X7 Balances for PGC
 PS R2 Precision Balance
 PS X2 Balances for PGC
 WLC Precision Balance
 PS X7 Precision Balance

AS 5Y Analytical Balance
 PS 5Y Balances for PGC
 PS R2.H Precision Balance
 WLC C/2 Precision balance
 PS R1 Precision Balance
 AS R2 PLUS Analytical Balance
 PS 5Y Precision Balance
 AS R1 PLUS Analytical Balance

Device dimensions

