



More information on the website  
[radwag.com/us/info,w1,W20](https://radwag.com/us/info,w1,W20)

# Under-hook weighing rack for AS and PS balances



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A rack for under-hook weighing intended for balances AS and PS.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field. The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

## Compatible with

AS X7 Analytical Balance  
 PS X2 Precision Balance  
 AS X2 PLUS Analytical Balances  
 PS X7 Balances for PGC  
 PS R2 Precision Balance  
 PS X2 Balances for PGC  
 WLC Precision Balance  
 PS X7 Precision Balance

AS 5Y Analytical Balance  
 PS 5Y Balances for PGC  
 PS R2.H Precision Balance  
 WLC C/2 Precision balance  
 PS R1 Precision Balance  
 AS R2 PLUS Analytical Balance  
 PS 5Y Precision Balance  
 AS R1 PLUS Analytical Balance

## Device dimensions

